

## Glycosyltransferase Affinity Gel - GDP

Description:	Guanosin-5'-diphosphat coupled with a spacer via pyrophosphate linkage to a carrier.
CDP-concentration:	> 10 $\mu$ mol GDP per ml gel slurry (method: UV-Spectroscopy after hydrolysis with 1 M HCl)
Application:	purification of fucosyltransferases by affinity chromatography
Order-no:	5085
Form of delivery:	grey-white suspension with 20% ethanol
Storage:	4°C, do not freeze!